


RoweSpike II



- For repeated withdrawal and aspiration of solutions, drugs and toxic substances
- Safety snap cap for simple one-hand handling and protection against contamination
- Easy handling due to ergonomic design
- Low priming volume
- Low piercing forces
- Easy admission
- Ventilation filter 0.1 µm PTFE membrane (hydrophobic bacteria filter)

Art.-No.	Material (body)	Filling volume (ml)	Connection	Hydrophilic filter	Cap	Specifics	PU ((pcs./box) / pal.)
A-6425	ABS - w	0.25	fLL	-	green	*	300 / 16,800
A-6426	ABS - w	0.25	fLL	5 µm	blue	*	300 / 16,800
A-6427	MABS - t	0.25	fLL	5 µm, chemo	red	*	150 / 16,800
A-6453	ABS - w	0.25	fLL	-	green	eco * / **	300 / 16,800
A-6454	ABS - w	0.25	fLL	5 µm	blue	eco * / **	300 / 16,800
A-6455	MABS - t	0.25	fLL	5 µm, chemo	red	eco * / **	150 / 16,800

w = white
t = transparent

* Safety Valve variants are available

** Deep-drawn fluid drain to empty the bottle completely



RoweMed AG - Medical 4 Life
www.rowemed.de • info@rowemed.de
Juri-Gagarin-Ring 4 • 19370 Parchim • Germany
Tel.: +49 (0) 3871 451280 • Fax : +49 (0) 3871 451282