

## RoweSpike II MK / Micro



- For repeated withdrawal of solutions, drugs and toxic substances and ventilation of vials and containers
- Easy handling due to ergonomic design

- Low priming volume
- Low piercing forces
- Ventilation filter 0.1 µm PTFE membrane (hydrophobic bacteria filter)

Art.-No.	Material (body)	Filling volume (ml)	Connection	Tip	Hydrophilic filter	Cap	PU ((pcs./box) / pal.)
A-6430	ABS - w	0.25	fLL	Metal cannula	5 µm	blue	300 / 16,800
A-6431	MABS - t	0.25	fLL	Metal cannula	5 µm, chemo	red	150 / 16,800
A-6433	MABS - t	0.25	fLL	Metal cannula	chemo, free flow	Safety Valve, green	150 / 16,800
A-6436	ABS - w	0.20	fLL	Micro / Metal cannula	5 µm	blue	300 / 16,800
A-6437	MABS - t	0.20	fLL	Micro / Metal cannula	5 µm, chemo	red	150 / 16,800
A-6444	ABS - w	0.25	fLL	Micro / Metal cannula	5 µm	Safety Valve, transparent	300 / 16,800
A-6449	MABS - t	0.25	fLL	Micro / Metal cannula	5 µm, chemo	Safety Valve, red	150 / 16,800

w = white

t = transparent

\* Safety Valve variants are available

RoweMed AG - Medical 4 Life  
 www.rowemed.de • info@rowemed.de  
 Juri-Gagarin-Ring 4 • 19370 Parchim • Germany  
 Tel.: +49 (0) 3871 451280 • Fax : +49 (0) 3871 451282

