

RowePaed M



- Integrated venting through a PTFE membrane
- Compatible with all Luer-Lock connectors
- Filter area 2.3 cm²
- Transparent housing
- bulk / non-sterile

Art.-No.	Material (body)	Filling volume (ml)	Connection	Hydrophilic Filter	Color Code	PU (pcs./box)
A-6510/04	PC	0.3	2x (mL / ID 2.0)	1.2 µm PET	yellow arrows	4,000
A-6508/04	PC	0.3	2x (mL / ID 2.0)	0.2 µm PES	-	4,000
A-6511/04	PC	0.3	2x (mL / ID 2.0)	0.2 µm nylon*	-	4,000
A-6512/04	PC	0.3	2x (mL / ID 2.0)	0.2 µm PET	-	4,000
A-6513/04	PC	0.3	2x (mL / ID 2.0)	1.2 µm PET	-	4,000
A-6514/04	PC	0.3	2x (mL / ID 2.0)	0.2 µm nylon*	red housing base	4,000
A-6515/04	PC	0.3	2x (mL / ID 2.0)	1.2 µm PET	blue housing base	4,000

* positively charged nylon membrane with endotoxin retention

RoweMed AG - Medical 4 Life
 www.rowemed.de • info@rowemed.de
 Juri-Gagarin-Ring 4 • 19370 Parchim • Germany
 Tel.: +49 (0) 3871 451280 • Fax : +49 (0) 3871 451282

