

## RoweBagSpike



- For repeated withdrawal and aspiration of solutions or drugs out of all conventional medical bags
- High security and maximum protection against contamination due to closed system
- Good priming characteristics

Art.-No.	Length (cm)	Tubing Ø (mm)	Material	Filling volume (ml)	Connection	Hydrophobic filter	Specifics	PU ((pcs./box) / pal.)
A-2500	6	3.0 x 4.1	PUR	0.9	fLL	1.2 µm	Pinch clamp	100 / 11,200
A-2501	6	3.0 x 4.1	PUR	0.9	fLL	1.2 µm	Safety Valve	100 / 11,200
A-2502	6	3.0 x 4.1	PVC	0.9	fLL	1.2 µm	Pinch clamp	100 / 11,200
A-2504	-	-	-	0.55	fLL	1.2 µm	Safety Valve	100 / 11,200
A-2505	4	3.0 x 4.1	PUR	0.3*	fLL	-	Safety Valve, drip chamber	50 / 5,600
A-2506	30	3.0 x 4.1	PVC	2.3	mLL	-	Safety Valve	50 / 5,600
A-2507	10	4.8 x 6.8	PVC	2.0	twist-off port	-	Pinch clamp	50 / 5,600

\* tube only

RoweMed AG - Medical 4 Life  
 www.rowemed.de • info@rowemed.de  
 Juri-Gagarin-Ring 4 • 19370 Parchim • Germany  
 Tel.: +49 (0) 3871 451280 • Fax : +49 (0) 3871 451282

