

## RowePhob 18 / 25



- Hydrophobic ventilation filter for safe retention of fluids, particles, bacteria and critically aerosols > 0.1 µm
- Good grip
- Good ventilation

- Compatible with all Luer-Lock connectors
- Different connectors

Art.-No.	Material (body)	Filling volume (ml)	Connection	Filter Ø (mm)	Filter	PU ((pcs./box) / pal.)
A-6250	MABS	0.2	fLL/mL	18	0.1 µm PTFE	1,000 / 16,000
A-6275	MABS	0.2	fLL/mLL	18	0.1 µm PTFE	1,000 / 16,000
A-6350	MABS	0.6	fLL/mL	25	0.1 µm PTFE	1,000 / 16,000
A-6375	MABS	0.6	fLL/mLL	25	0.1 µm PTFE	1,000 / 16,000

RoweMed AG - Medical 4 Life  
 www.rowemed.de • info@rowemed.de  
 Juri-Gagarin-Ring 4 • 19370 Parchim • Germany  
 Tel.: +49 (0) 3871 451280 • Fax : +49 (0) 3871 451282

