

## RowePaed



- Especially for paediatrics
- Integrated venting through a PTFE membrane
- Compatible with all Luer-Lock connectors
- Can be connected to additional systems easily
- Visual drug flow control
- Rotating lock nut
- Filter area 4 cm<sup>2</sup>

Art.-No.	Material (body)	Filling volume (ml)	Connection	Filter	PU ((pcs./box) / pal.)
A-6600	MABS	0.7	fLL/mLL	0.2 µm nylon*	1,000 / 16,000
A-6601	MABS	0.7	fLL/mLL	1.2 µm PET	1,000 / 16,000
A-6602	MABS	0.7	fLL/mLL	0.2 µm PET	1,000 / 16,000
A-6603	MABS	0.7	fLL/mLL	0.2 µm PET	1,000 / 16,000
A-6604	MABS	0.7	fLL/mLL	1.2 µm PET	1,000 / 16,000
A-6605	MABS	0.7	fLL/mLL	0.2 µm nylon*	1,000 / 16,000
A-6606	MABS	0.7	fLL/mLL	0.2 µm PES	1,000 / 16,000

\* positively charged nylon membrane with endotoxin retention

♦ without closing elements

RoweMed AG - Medical 4 Life  
 www.rowemed.de • info@rowemed.de  
 Juri-Gagarin-Ring 4 • 19370 Parchim • Germany  
 Tel.: +49 (0) 3871 451280 • Fax : +49 (0) 3871 451282

