


RoweInfuset with RowePaed M



- Infusion set with paediatric filter and transparent tubings
- Integrated venting through a PTFE membrane
- Compatible with all Luer-Lock connectors
- Can be connected to additional systems easily
- Filter area 2.3 cm²
- Transparent connectors / tube / housing
- Visual drug flow control
- PC housing

Art.-No.	Material (tube)	Length (cm)	Tubing Ø (mm)	Filling volume (ml)	Connection	Filter	PU ((pcs./box) / pal.)
A-2310	PUR	2 x 3	1.0 x 2.0	0.5	fLL/mL + rotating lock nut	0.2 µm nylon*	800 / 12,800
A-2369	PUR	2 x 10	1.0 x 2.0	0.6	fLL/mLL	0.2 µm nylon* + slide clamp	800 / 12,800
A-2370	PUR	3 x 5	1.0 x 2.0	0.9	fLL/mLL+ rotating lock nut	0.2 µm nylon* + check valve	800 / 12,800
A-2372	PUR	2 x 10	1.0 x 2.0	0.6	fLL/mLL	1.2 µm PET + slide clamp	800 / 12,800
A-2376	PUR	2 x 10	1.0 x 2.0	0.6	fLL / mLL	1.2 µm PET + slide clamp + color code: blue	800 / 12,800

* positively charged nylon membrane with endotoxin retention

RoweMed AG - Medical 4 Life
 www.rowemed.de • info@rowemed.de
 Juri-Gagarin-Ring 4 • 19370 Parchim • Germany
 Tel.: +49 (0) 3871 451280 • Fax : +49 (0) 3871 451282

