

RowePhil 3 / 18



RowePhil 3:

- Filter Straw for withdrawal and filtration of solutions and drugs from glass ampoules
- Compatible with all Luer-Lock connectors
- Visual flow control
- Needle free withdrawal
- Flexible plastic straw
- Customized systems are available
- MABS housing
- Filter diameter 3 mm



RowePhil 18 (Cannula / Spike):

- Hydrophilic injection filter for protection against particle contamination of small amounts of drugs
- Low priming volume

Art.-No.	Length (cm)	Tubing Ø (mm)	Filling volume (ml)	Connection	Filter	Material	Specifics	PU ((pcs./box) / pal.)
A-1312	10	1.5 x 2.7	~ 0.11	fLL	-	MABS / PUR	Withdrawal Straw	1,000 / 16,000
A-2390 *	10	1.5 x 2.7	0.25	fLL	5 µm PET	MABS / PUR	Filter Straw	1,000 / 16,000
A-2391 *	4.5	1.5 x 2.7	0.13	fLL	5 µm PET	MABS / PUR	Filter Straw	1,000 / 16,000
A-2392 *	-	-	-	fLL/mLL	5 µm PET	-	Filter	on request
A-6230	-	0.2	-	fLL/Spike	5.0 µm PET (hydrophil)	MABS	-	1,000 / 16,000
A-6276 **	-	0.2	-	fLL/Cannula 18 G	0.1 µm PTFE (hydrophob)	MABS	-	1,000 / 16,000
A-6284	-	0.2	-	fLL/Cannula 18 G	5.0 µm PET (hydrophil)	MABS	-	1,000 / 16,000

* reduces risk of glass particle contamination
 ** hydrophobic ventillation filter

RoweMed AG - Medical 4 Life
 www.rowemed.de • info@rowemed.de
 Juri-Gagarin-Ring 4 • 19370 Parchim • Germany
 Tel.: +49 (0) 3871 451280 • Fax : +49 (0) 3871 451282

