

## RoweSpike I



- For repeated withdrawal of solutions, drugs and toxic substances and ventilation of vials and containers
- Safety snap cap for simple one-hand handling and protection against contamination
- Easy handling due to ergonomic design
- Low priming volume
- Low piercing forces
- Ventilation filter 0.1 µm PTFE membrane (hydrophobic bacteria filter)

Art.-No.	Material (body)	Filling volume (ml)	Connection	Hydrophilic filter	Cap	Specifics	PU ((pcs./box) / pal.)
A-6400	ABS - w	0.25	fLL	-	green	*	300 / 16,800
A-6401	ABS - w	0.25	fLL	5 µm	blue	*	300 / 16,800
A-6402	MABS - t	0.25	fLL	5 µm, chemo	red	*	150 / 16,800
A-6407	ABS - w	0.22	fLL	-	green	eco * / **	300 / 16,800
A-6408	ABS - w	0.22	fLL	5 µm	blue	eco * / **	300 / 16,800
A-6409	MABS - t	0.22	fLL	5 µm, chemo	red	eco * / **	150 / 16,800

w = white  
t = transparent

\* Safety Valve variants are available

\*\* Deep-drawn fluid drain to empty the bottle completely



RoweMed AG - Medical 4 Life  
www.rowemed.de • info@rowemed.de  
Juri-Gagarin-Ring 4 • 19370 Parchim • Germany  
Tel.: +49 (0) 3871 451280 • Fax : +49 (0) 3871 451282

