

RoweSpike VI SnapFit



- For repeated withdrawal of solutions, drugs and toxic substances and ventilation of vials and containers
- Safety snap cap for simple one-hand handling
- Easy handling due to ergonomic design
- Protection against contamination through safety valve (SV)



- Low piercing forces
- Ventilation filter 0.1 µm PTFE membrane (hydrophobic bacteria filter)
- Safe vial connection for 20 mm bottle
- Deep-drawn fluid drain to empty the bottle completely

Art.-No.	Material (body)	Filling volume (ml)	Connection	Tip	Hydrophilic filter	Specifics	PU ((pcs./box) / pal.)
A-6472	MABS	0.13	fLL	Spike	-	-	150 / 16,800
A-6473	MABS	0.20	fLL	Spike	-	Safety Valve	150 / 16,800
A-6485	MABS	0.25	fLL	Spike	5 µm, chemo	red cap	150 / 16,800
A-6486	MABS	0.25	fLL	Spike	5 µm, chemo	Safety Valve, red	150 / 16,800

RoweMed AG - Medical 4 Life
www.rowemed.de • info@rowemed.de
 Juri-Gagarin-Ring 4 • 19370 Parchim • Germany
 Tel.: +49 (0) 3871 451280 • Fax : +49 (0) 3871 451282

