


RoweSpiralLines / RoweTransfer-Sets



- pressure resistance up to 21 bar (305 psi)
- standard luer connectors

- visual drug flow

Art.-No.	Material (tube)	Length (cm)	Tubing Ø (mm)	Specifics	Connection	PU ((pcs./box) / pal.)
A-7003	PVC	150	1.65 x 3.15	-	fLL/mLL	50 / 3,600
A-7009	PVC	180	1.65 x 3.15	-	fLL/mLL	50 / 3,600
A-7012	PVC	-	3.00 x 4.10 2.50 x 4.10	RoweTransfer-Set Tip: drip chamber 2 x Y-connector 2 x check valve	fLL/mLL	50 / 3,600
A-7021	PVC	180	1.65 x 3.15	check valve, Y-connector	fLL/mLL	50 / 3,600
A-7031	PVC	-	3.00 x 4.10 2.50 x 4.10	RoweTransfer-Set Tip: Spike check valve	fLL/mLL	50 / 3,600
A-7039	PVC	150	1.65 x 3.15	check valve, Y-connector	fLL/mLL	50 / 3,600
A-7040	PVC	210	1.65 x 3.15	check valve, Y-connector	fLL/mLL	50 / 3,600
A-7041	PVC	210	1.65 x 3.15	-	fLL/mLL	50 / 3,600

customized systems are available

RoweMed AG - Medical 4 Life
 www.rowemed.de • info@rowemed.de
 Juri-Gagarin-Ring 4 • 19370 Parchim • Germany
 Tel.: +49 (0) 3871 451280 • Fax : +49 (0) 3871 451282