

## RowePhil 3 / 18



### RowePhil 3:

- Filter Straw for withdrawal and filtration of solutions and drugs from glass ampoules
- Compatible with all Luer-Lock connectors
- Visual flow control
- Needle free withdrawal
- Flexible plastic straw
- Customized systems are available
- MABS housing
- Filter diameter 3 mm



### RowePhil 18 (Cannula / Spike):

- Hydrophilic injection filter for protection against particle contamination of small amounts of drugs
- Low priming volume

| Art.-No.  | Length (cm) | Tubing Ø (mm) | Filling volume (ml) | Connection       | Filter                  | Material   | Specifics        | PU ((pcs./box) / pal.) |
|-----------|-------------|---------------|---------------------|------------------|-------------------------|------------|------------------|------------------------|
| A-1312    | 10          | 1.5 x 2.7     | ~ 0.11              | fLL              | -                       | MABS / PUR | Withdrawal Straw | 1,000 / 16,000         |
| A-2390 *  | 10          | 1.5 x 2.7     | 0.25                | fLL              | 5 µm PET                | MABS / PUR | Filter Straw     | 1,000 / 16,000         |
| A-2391 *  | 4.5         | 1.5 x 2.7     | 0.13                | fLL              | 5 µm PET                | MABS / PUR | Filter Straw     | 1,000 / 16,000         |
| A-2392 *  | -           | -             | -                   | fLL/mLL          | 5 µm PET                | -          | Filter           | on request             |
| A-6230    | -           | 0.2           | -                   | fLL/Spike        | 5.0 µm PET (hydrophil)  | MABS       | -                | 1,000 / 16,000         |
| A-6276 ** | -           | 0.2           | -                   | fLL/Cannula 18 G | 0.1 µm PTFE (hydrophob) | MABS       | -                | 1,000 / 16,000         |
| A-6284    | -           | 0.2           | -                   | fLL/Cannula 18 G | 5.0 µm PET (hydrophil)  | MABS       | -                | 1,000 / 16,000         |

\* reduces risk of glass particle contamination  
 \*\* hydrophobic ventillation filter

RoweMed AG - Medical 4 Life  
 www.rowemed.de • info@rowemed.de  
 Juri-Gagarin-Ring 4 • 19370 Parchim • Germany  
 Tel.: +49 (0) 3871 451280 • Fax : +49 (0) 3871 451282